



SPECIALS

OCT / NOV / DEC 2018

COMPOSITES | GLASS IONOMERS
TOOTH WHITENING | EQUIPMENT
CEMENTS | ADHESIVES | RIVA STAR
ETCHANTS | SEALANTS
AMALGAMS | ACCESSORIES

WHITENING

ALL POLA KITS

3+1*

BUY 3, GET 1 FREE

5+2*

BUY 5, GET 2 FREE

10+7*

BUY 10, GET 7 FREE



POLA DAY

HYDROGEN PEROXIDE BASED ADVANCED TRAY TOOTH WHITENING SYSTEM

- High water content
- From 30 minutes once a day
- Contains fluoride
- Available in 3% and 6% hydrogen peroxide
- Neutral pH and desensitising additives
- Stays in place
- Pleasant spearmint flavoured gel

POLA NIGHT

CARBAMIDE PEROXIDE BASED ADVANCED TRAY TOOTH WHITENING SYSTEM

- High water content
- From 90 minutes once a day
- Contains fluoride
- Available in 10% and 16% carbamide peroxide
- Neutral pH and desensitising additives
- Stays in place
- Pleasant spearmint flavoured gel

ALLOYS

GS80 (ADMIX & SPHERICAL), PERMITE & ULTRACAPS

3+1*

BUY 3, GET 1 FREE ON 500 CAPS

5+1*

BUY 5, GET 1 FREE ON 50 CAPS



GS-80

ECONOMICAL NON-GAMMA 2 ADMIX ALLOY

- High early compressive strength
- Great corrosion resistance
- Excellent marginal sealing
- Optimal handling characteristics
- Dimensional stability
- Over 400 Million restorations worldwide



PERMITE

HIGH SILVER NON-GAMMA 2 ADMIX ALLOY

- Superior handling characteristics
- Void-free restorations with no marginal leakage.
- Extremely high resistance to corrosion
- Mechanical and chemical stability
- High strength and hardness for extreme wear resistance



COMPOSITES



The Dental Advisor, Vol. 33,
No. 07, September 2016

AURA BULK FILL

THE SIMPLICITY OF BULK FILLING

- One Bulk Fill placement up to 5mm depth of cure†
- 20 sec. LED light cure
- Natural-looking finish with just one universal shade
- High depth of cure without the need for a capping layer
- Easy placement in a matter of seconds
- The waxy consistency is packable and non-sticky to optimise handling

1+1* BUY 1, GET 1 FREE

ETCHANT



SUPER ETCH

37% PHOSPHORIC ACID ETCH GEL

- Stays where placed.
- Does not run onto soft tissue
- Washes away easily with water
- Easily controlled gel which flows under direction, not gravity
- Blue tint for easy placement visibility
- Thixotropic gel to improve handling

2+1* BUY 2, GET 1 FREE

GLASS IONOMERS



RIVA SELF CURE

REGULAR & HIGH VISCOSITY, SELF-CURING, TOOTH
REMINERALISING GLASS IONOMER RESTORATIVE MATERIAL

- Very high fluoride releasing and recharging ability to create more acid resistant restored teeth



RIVA LIGHT CURE

REGULAR & HIGH VISCOSITY, LIGHT
CURED, RESIN REINFORCED GLASS
IONOMER RESTORATIVE MATERIAL

- Orange light safe tip prevents premature setting

**BUY 7 BOXES
OF RIVA
GET ULTRAMAT S
FOR FREE**



LUTING CEMENT



RIVA CEM

SELF-CURING RESIN MODIFIED GLASS IONOMER
LUTING CEMENT

- Easy to mix paste/paste system
- Light tack cure option for easy removal of excess paste
- High bond strength
- BPA free
- Greater stain resistance
- No wastage at tip due to hand mix

1+1* **RIVA CEM**
**BUY 1,
GET 1 FREE**

PARTICIPATING DEALERS

BDSI	0845 094 2722	Survival	32 0118 951 6161
BF Mulholland	0289 445 2668	The Dental Directory	0800 585 586
Budget	08700 100 770	Trycare Dental	0800 269 962
CTS Dental Supplies	0173 776 5400	Wright Cottrell	0800 668 899
Dental Health Boutique	08456 017 086		
Dental Sky	0800 294 4700		
Henry Schein	08700 102043		
Kent Express	01634 878787		

*Valid on invoices dated 1st October - 31st December, 2018. Please allow 2-4 weeks for delivery. Offers cannot be combined. All claims for free goods must be made within 60 days from the end of promotion date. Free goods will be of equal or lesser value. Offers are valid in the UK only.

SDI Germany GmbH

Email uk@sdi.com.au or
contact **Kyra Hillebrand** on
08000 022 55 734
for details of your nearest
SDI representative

WWW.SDI.COM.AU
WWW.POLAWHITE.COM.AU



YOUR
SMILE.
OUR
VISION.